



Biological (Far-field) Microscopy Beyond the Diffraction Limit of Light

Joerg Bewersdorf

Departments of Cell Biology and of Biomedical Engineering



1



2


The Diffraction Limit



Ernst Abbe
(1872)

diffraction limit

structures smaller than
half a wavelength
cannot be resolved.



cf. Airy disk radius:

$$r_0 \approx 0.61 \frac{\lambda}{n \sin \alpha}$$

wavelength
 λ

refractive index
 n
aperture angle
 α

numerical aperture (NA)

3

Resolution in Far-field Microscopy

$$r_0 \approx 0.61 \frac{\lambda}{n \sin \alpha}$$


	λ	d_{\min}
light	~ 500 nm	~ 250 nm
X-ray	~ 2 nm	~ 25 nm
electron	~ 0.001 nm	~ 0.1 nm (>2 nm)

4

Why Fluorescence Light Microscopy?

Advantages



- Wide variety of labels and targets
 - Specific
 - Multi-color
 - Sensitive
- Imaging live specimens
- 3D imaging inside samples
- Easy sample preparation compared to EM



5


A Fascinating World Below 100 nm

cristae
(Essential Cell Biol. 2010)






→ Need for better resolution in light microscopy


synaptic vesicle
(Takamori et al., Cell 2006)



phage
(Bebeacua et al.; J Virol 2012)




6



7

Overview of Super-resolution (Far-field Fluorescence) Microscopy

Light	<p>Pushing the resolution within the diffraction limit</p> <ul style="list-style-type: none"> • 4Pi microscopy • Structured Illumination Microscopy (SIM) • Image-scanning microscopy, iSIM, AiryScan, ... 	<p>Resolution improvement</p> <p>} ~2x</p>
Sample-light interaction	<p>Breaking the diffraction limit of light</p> <ul style="list-style-type: none"> • Spatially targeted off-switching techniques (STED, ...) • Single-molecule localization techniques (PALM/STORM/PAINT/... and MINIFLUX) 	<p>} ~5-10x (50x)</p> <p> 2014</p>
Sample	<p>Making the sample bigger</p> <ul style="list-style-type: none"> • Expansion Microscopy 	<p>} ~4-20x</p>


8

STED Microscopy (or: STED 'Nanoscopy')

Breaking the diffraction limit by
targeted off-switching of molecules


9

Principle of STED Nanoscopy




See Hell & Wichmann, *Opt Lett* 1994

10




11



12


A Custom-built STED Nanoscope for Live-cell Imaging



inspired by: Laser lines: Gottfert *et al.*, Biophysical Journal (2013)
 SLM configuration: M.O Lenz *et al.*, Journal of Biophysics. (2014)


13

Example of STED Microscopy




Travis Gould

14




15



16

Dynamics of Micro-fenestrations in ER Sheets



Lena Schroeder

COS7 cell at 37 °C

SNAP-KDEL BG-SiR

4 frames/second

Bleaching compensation
applied

Collaboration with Shirin
Bahmanyar

17

3D STED Microscopy for Deep-tissue Live-animal Imaging


Challenge: Scattering and wavefront aberrations

- Solution:**
- 2-Photon Excitation**
→ reduced scattering (proportional to λ^{-3} to λ^{-4})
 - Far-red Dyes**
→ reduced scattering in detection beam path
 - Adaptive Optics**
→ reduced optical aberrations
 - Long Working Distance Objective**

Velasco et al., Optica 2021

18

2-photon STED Microscopy




Moneron & Hell, *Opt. Express* (2009)
 Ter Veer, ..., Nagerl, *Methods Mol. Biol.* (2017)
 Velasco et al., *Opt. Lett.* (2015)

Mary Grace Velasco

19

Aberration Correction




Review: Booth, *Light: Science & Applications* (2014)

Mary Grace Velasco

20

Aberration Correction




Aviles-Espinosa et al., *Biomed. Opt. Express* (2011)

Wang et al., *Nature Methods* (2014)

Mary Grace Velasco

21

Aberration Correction




Mary Grace Velasco

22


Custom-built 2-photon STED Microscope

①
Spatial Light Modulator (SLM)



... for 3D Depletion


*Lenz et al., *J. Biophotonics* (2014)



Velasco et al., *Optica* 2021


23

Custom-built 2-photon STED Microscope




Velasco et al., *Optica* 2021

24



25




26


3D STED Microscopy In Mouse Brain Tissue Section

Combination: STED Microscopy
2-Photon
Adaptive Optics

2-Photon



2-Photon STED



Astrocytes in mouse brain section
Anti-GFAP ATTO594


Unpublished
Mary Grace Velasco, Mark Lessard
Bewersdorf Lab

27

Single-molecule Localization Microscopy (PALM/STORM/FPALM/...)


Breaking the diffraction limit by
stochastically switching of individual molecules

28



29

Single-Molecule Localization



For isolated molecules, determine the center of the diffraction-limited image to find its position.

Typically done by fitting a model function to the image.

The precision depends primarily on S/N .


Simple approximation:
$$\sigma_{loc} \approx \frac{\sigma_{PSF}}{\sqrt{N}}$$

More accurately:

$$\sigma_{loc} \approx \sqrt{\frac{\sigma_{PSF}^2 + a^2/12}{N} \left(1 + 4\tau + \sqrt{\frac{2\tau}{1+4\tau}} \right)}$$


with:
$$\tau = \frac{2\pi b(\sigma_{PSF}^2 + a^2/12)}{Na^2}$$

σ_{PSF} : PSF width
 N : signal photons
 a : pixel size
 b : background photons per pixel



Rieger & Stallinga ChemPhysChem 2014

30



31

3D Single-molecule Switching Nanoscopy

3D SMS resolution:

XY	~20-30 nm
Z	~50-80 nm

4Pi-SMS/iPALM:

XY	~20 nm
Z	~10-20 nm

iPALM: Shtengel *et al. PNAS*, 2009
4Pi-SMS: Aquino *et al. Nature Methods*, 2011

~~Reason: elongated PSF~~

Solution: 4Pi detection

S.W. Hell 1992
Gugel, Bewersdorf *et al. Biophys. J.*, 2004

Astigmatism:
Huang *et al. Science*, 2008
Biplane:
Juette *et al. Nature Methods*, 2008
Double Helix:
Pavani *et al. PNAS*, 2009
Phase Ramp:
Baddeley *et al. Nano Research*, 2011
...

32

4Pi-SMS Nanoscopy

Interferometric Images

$\phi=0$
 $\phi=\pi/2$
 $\phi=\pi$
 $\phi=3\pi/2$

Intensity
Fluorophore Z position

George Sirinakis

Shtengel, G. et al. *PNAS* **2009**, 9, 3125-30
Aquino, D. et al. *Nat. Methods* **2011**, 8, 353-359

33

4Pi-SMS Nanoscopy


20 x 20 x 10 nm resolution
Adaptive optics enables >5 μm thick samples

George Sirinakis, Edward Allgeyer

34

4Pi-SMS Nanoscopy


20 x 20 x 10 nm resolution
Adaptive optics enables >5 μm thick samples



George Sirinakis, Edward Allgeyer


35

4Pi-SMS raw data



36

The Golgi Complex Resolved by Light Microscopy



GRASP65

HeLa cell

Collaboration with
Rothman Lab

Yongdeng Zhang, Lena Schroeder et al., Nature Methods 2020

37


Expansion Microscopy

Enlarge the sample to resolve
originally small structures

38

Expansion Microscopy


- Invented 2015
 - Chen, Tillberg, Boyden, Science 347:543-548
- Sample is physically expanded 4-5x in a hydrogel
- Initially: only specific proteins
 - Probes cross-linked to gel
 - all proteins degraded
 - requires labeling before expansion process, which is not ideal (crowded environment; labels also expand)
- 2016: protease-free approach (MAP, U-ExM)
 - Proteins not degraded anymore
 - post-expansion labeling possible
- 2017: iterative expansion (i-ExM)
 - Up to 20x (~4.5 x 4.5) expansion
 - makes very small structures visible; but degraded all proteins (see problems above)




Wassie, Zhao & Boyden, Nat. Meth. 2019

39

Expansion Microscopy Examples




mouse hippocampus expressing virally delivered Brainbow3.0 epitopes




smFISH of NEAT1 lncRNA in the nucleus of a HeLa cell

Wassie, Zhao & Boyden, Nat. Meth. 2019


40




41



42




43



44

Pan-ExM reveals centriole ultrastructure




63X/1.2 NA water objective; confocal microscopy

M'Saad & Bewersdorf, *Nat. Comm.* (2020)

45

Synaptic Organization in Cultured Neurons



M'Saad et al., *unpublished*

46

Acknowledgements

Current:

Kenny Chung
Lukas Fuentes
Kevin Hu
Phylicia Kidd
Xiangzhou Lao
Dong-Ryoung Lee
Yang Li
Zach Marin
Hannahmariam Mekbib
Ons M'Saad
Florian Schueder
Yuan Tian

Former:

Andrew Barentine
Edward Courchaine
Xiang Hao
Fang Huang
Emil Kromann
Mark Lessard
Yu Lin
Timy Phan
Bennett Rollins
Lena Schroeder
Helen Sun
Mary Grace Velasco
Yongdeng Zhang

Cell Biology, Yale

Jim Rothman
Felix Rivera-Molina
Derek Toomre
Jeeyun Chung
Pietro de Camilli
Zhao Zhang
Nathan Williams
Yang Yang
Chenxiang Lin

Mol., Cell. and Dev. Biology, Yale

Tom Pollard
Holly Merta
Shirin Bahmanyar

Microbial Pathogenesis, Yale

Maria Lara-Tejero
Jorge Galan

Genetics, Yale

Miao Liu
Siyuan (Steven) Wang
Mark Pownall
Antonio Giraldez

Mol. Biophysics & Biochem., Yale

Tony Koleske
Karla Neugebauer

Neurology, Yale

Ilona Kondratiuk
Thomas Biederer

University of Auckland, NZ

David Baddeley

Cambridge University

Daniel St Johnston
Edward Allgeyer
George Sirinakis

Oxford University

Martin Booth
Jingyu Wang
Jacopo Antonello
Jens Rittscher
Ian Dobbie

UT Southwestern

Kevin Dean
Reto Fiolka
Bo-Jui Chang

Janelia Research Campus

Jonathan Grimm
Luke Lavis

EMBL

Yiming Li
Jonas Ries

Supported by
wellcome trust



National Institutes
of Health



Disclaimer: J.B. has significant financial interest in Bruker Corporation and Hamamatsu Photonics and a startup company in the Expansion Microscopy field.